DE

EN

b-intense®

Bedienungsanleitung (instruction for use, IFU) b-intense® MED

inkl. Anschlussplan b-control MED
(incl. connection plan b-control MED)

INHALT

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CONTENT

Page	Topic
33	General safety guidelines
34	General safety notices
38	General warning notices and recommendations
40	b-intense ® MED: transportation and storage
41	General Application Notes
42	Scope of application of b-intense ® MED
45	Information on radiation emissions
46	b-control MED: Operation
48	b-control MED: Service menu
49	Cleaning and care of the b-intense ® MED
49	General cleaning of the b-intense ® MED
50	Maintenance instructions
50	Disposal of the b-intense ® MED
51	Certificate of Warranty
52	b-intense MED: marking plates
55	b-intense MED: Locations of EMISSION APERTURES
56	b-intense MED: technical data
58	Connection plan b-control MED control unit
59	Appendix A: Informations reg. EMF according to EN 60601-1-2
61	Appendix B: List of replaceable cables

IMPORTANT SAFETY GUIDELINES. PLEASE READ CAREFULLY AND KEEP SAFE FOR LATER USE.

General safety guidelines

All electrical components comply with safety regulations EN60601, EN62471 and EN62304. Please check the **b-intense**[®] infrared heaters before use for any foreign bodies to prevent accidents. The **b-intense**[®] MED complies with the relevant safety guidelines for medical devices. In case of any kind of damage to the heaters, do not start the infrared cabin and directly contact **b-intense**[®].

b-intense[®]

Manufacturer: **b-intense**[®] gmbh, branch office of CB-Capital Beteiligungs GmbH
Am Jungferenberg 17, A-2201 Gerasdorf by Vienna, Austria
office@b-intense.at

CE 1304

General warning sign

Follow the instruction for use

Protective earth (earth)

Warning, electricity

Caution, radiation emitting device

Aperture for optical radiation

Hot surface. Do not touch the grid of the front and back heater.

Pushing prohibited. Do not move the **b-intense**[®] MED.

General safety notices

Observe the safety instructions before use.

Store and operate in a dry environment.
Humidity less than 60% for operation.
Not condensing.

Safety Warning: Never cover an infrared heater!
Make sure that the radiator is not covered and no items are put on it.
Covering the infrared radiator may cause a fire hazard.

Follow all safety instructions in this manual!
Work on the electrical wiring must only be carried out by a qualified electrician. Turn the b-intense® MED 110V off first and always unplug it before repair and maintenance or cleaning work is carried out!
There is a danger of electric shock!

These operating instructions contain specific information about the radiation emissions, about options for radiation protection for users and about false use and risks through the installation (see chapter "Information on radiation emission").

IR emitted from this device may cause eye injury
Avoid eye exposure
Note: Skin must be free during treatment and need not to be protected!

Aperture for optical radiation

Device excess Risk Group 1 at distance 400 mm or less

Check the required power supply. The operating voltage is announced on the marking plate of the b-intense MED (see page 80).

To avoid the risk of electric shock, this device must only be connected to a mains supply with a protective earth conductor.

The accessibility of the mains plug must always be ensured for disconnecting from the mains supply in an emergency.

Only original b-intense® component parts are safe for use. The use of not certified components can lead to damages or injuries. Warranties and guarantees expire with the use of not certified components.

1. The used socket should be secured with a fault-current circuit breaker (FI) with a trip current of max 30mA and a circuit breaker (safety device) of max 16A.
2. The electrical installation must be checked annually by a licensed electrical company.
3. The IR cabin must not be operated on a multiple socket (fire hazard).
4. The used socket should be easily accessible after installation.
5. The IR cabin must only be operated on an earthed shockproof socket, which was installed and checked by a licensed electrical company.
6. Only the original power cable with the length of 3 m must be used for connection to the power supply.
7. The power cable is not allowed to be extended.

Plugging the cabin in and out is only permitted by pulling the plug, which is connected to the cable. Don't pull the cable itself, this could damage the cable and subsequently lead to an electric shock.

Using this device next to other devices should be avoided, for it could lead to incorrect operation. If it is necessary, all affected devices should be monitored to secure regular operation.

To connect the electrical connection lines correctly, use the connection plan.

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the b-intense® MED, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

If the Essential Performance is lost or degraded due to EM Disturbances, the b-intense® MED will switch off all infrared heaters and the b-control MED unit will switch to the Stand-by mode. The treatment cannot be carried out.

The b-intense® MED is only intended to be used in professional healthcare environment!

Visually check the b-intense® MED regularly. In case of visible damage unplug the b-intense® MED from the power supply and contact us via office@b-intense.at.

Before working on the controller, the power supply must be switched off on all poles (pull the main plug) and prevented it from being switched on.

The use of components not licensed by b-intense® can cause higher electromagnetic emissions or reduce the electromagnetic immunity of the device and lead to improper operation and is therefore prohibited.

During regular operation, the IR heaters can reach high surface temperatures, which can cause burns when touched.

Only original b-intense® components (IR heaters, fan, sensor, etc.) must be connected to the control system. Spare parts have to be ordered from b-intense®.

Visually check the cabins installation before the first use (there should be no brittle wood, no broken glass, no nails or screws sticking out, the cabin should stand stable and not loose).

The safe working load for seat and cabin is 120 kg.

Modification of the b-intense® MED is not permitted.

Hot surface. Do not touch the grid of the front and back heater.

Pushing prohibited. Do not move the b-intense® MED.

If any malfunctions of the firmware are noticed, please unplug the b-intense® MED and report immediately to office@b-intense.

General warning notices and recommendations

Make sure that liquids do not come into contact with the heaters or with the timber facings.

Ensure that the heater is not covered and that there are no objects placed on it.

Avoid any contact with the heater while it is in operation and immediately afterwards. The high surface temperature of the heater can cause burns.

Ensure that your hair does not touch the heater while it is in operation. Tie your hair up to have your back completely free.

Allow the infrared heater to cool down for 9 minutes before re-starting the cabin.

Children are not allowed to play with the device or to carry out any repairs, cleaning or maintenance work.

Do not stare at an infrared heater for more than 9 minutes. We recommend safety glasses for sensible eyes (at least protection level 4-4 by EN 171).

Users with limited sensation of heat and users under the effects of alcohol, drugs or tranquilizers are not allowed to use the cabin.

Do not use the infrared cabin within 24 hours after UV irradiation from an artificial source (e. g. solarium) or sunbathing.

Reduce the intensity of the back heater immediately, if you feel burning or pain in the radiated area.

We recommend a warm shower and a nurturing, moisturizing lotion after the infrared application.

If you experience any discomfort, stop your infrared treatment immediately! If overheating occurs or if you feel the temperature is too high, you can reduce the radiation intensity or open the door and switch on the ventilation.

Only exchange electrical components with original components from the manufacturer. All repairs must be performed by an engineer qualified by **b-intensive®**. The electrical components must not be opened.

Consult your doctor on your own responsibility concerning your infrared treatment after operations, in case of illness, injuries, pregnancy or in the case of implants (e. g. after plastic or orthopedic surgery), also in any case that involves a risk of overheating, such as users with cardiovascular diseases or while taking medication.

Infrared used by children

Children may not use the infrared cabin.

Special note: Red back

A distinct redness in the exposed back area is normal and desired, and happens due to the increased blood circulation. If the redness is still visible more than 3 hours after application, reduce the radiation intensity the next time. We advise to contact a doctor, should the redness still be visible after 12 hours.

Warning

In the event of malfunctions, damages or injuries during application, please report immediately to **office@b-intense.at**.

Interruption of treatment

If the treatment is interrupted due to a defect, it should be restarted immediately in order to fulfil the treatment indications.

b-intense® MED: transportation and storage

Transportation of b-intense® MED

The **b-intense®** MED must be transported on pallets by a transport company suitable for pallet transport and fragile goods. In case of visible damage to the packaging, please contact the **b-intense®** office immediately. The cabin must not be assembled until release.

Storage of the b-intense® MED

The storage of the **b-intense®** MED should be as short as possible and not exceed one month in order to be able to complete the **b-intense®** quality check quickly. The stored product must remain on the pallet until the final installation.

Transport and storage conditions

Ambient Temperature: 0-40 °C

Relative humidity: 10-100 %

Atmospheric pressure: 50-110 kPa

Placement of the b-intense® MED

The best location is close to a socket, which will also be easily accessible after the assembly of the cabin. Ensure that there are no tripping hazards near the **b-intense®** MED.

Please consider the following recommendations for installation and storage:

- › The cabin must be stored and installed indoors.
- › The ambient temperature of the **b-intense®** MED should be at least 18 °C (room temperature) and at maximum 35 °C.
- › The optimal humidity should not exceed 60 %. If this is not possible, suitable ventilation of the space to be used must be arranged.
- › The cabin must be protected from any direct water spray and from wet or damp from the floor. Condensation in the room where the cabin is sited must be prevented and adequate ventilation must be ensured. Condensation on the components of the cabin must be avoided.
- › Ensure that the floor surface for the cabin is flat and level. Minor irregularities in the surface can be offset by the adjustable feet fixed to the underside of the cabin base.
- › The space between the wall and the cabin of at least 2 cm has to be considered.
- › The distance from the cabin roof to the ceiling above should be at least 10 cm.

The accessibility of the mains plug must always be ensured for disconnecting from the mains network in an emergency.

Please observe the following safety-relevant parameters:

- › The horizontal minimum distance of the back-heater to burning materials is 10 cm.
- › The horizontal minimum distance of the front-heater to burning materials is 30 cm.
- › The heaters must not be covered.
- › The power circuit must be secured using a fault-current circuit breaker (FI) with a trip current of 30 mA.

Service and Maintenance

Annual preventive replacement (or replacement after 100 hours of operation) of two 10 A fuses inside of the main unit by service personnel is recommended.

Service and maintenance has to be performed by service personnel. Users must not do any service and maintenance. The expected life time of the system is 10 years.

General Application Notes

Visually check the **b-intense**® MED for any damage before every use.

Please only enter the cabin with bathing shoes for hygienic reasons.

We recommend infrared sessions of approx. 30–50 minutes, depending on medical advice. After 60 minutes the cabin switches off automatically.

A **b-intense**® MED needs no warming up time. We still recommend warming the cabin up to 28 – 2 °C interior temperature for optimal results.

Mind your fluid balance! Drink a glass of water before and after each use of the infrared cabin.

The temperature in the **b-intense**® MED reaches between 30 °C and 40 °C.

The use of a towel while sitting in the cabin is recommended for hygienic reasons.

After use, a lukewarm shower is recommended.

Disinfect the seat, back and head rest according to prescription (p. 50) after every infrared application.

The infrared radiation must strike the skin directly. The application must occur with the upper body naked and leaned back, so that your spine is positioned directly in front of the back radiator. Adjust the neck support prior to application, so that the radiation does not hit the back of your head directly. Therefore grab the neck support with both hands central on the right and on the left side and push it into the desired position slowly with steady force.

Regulate the settings of the **b-intense**® MED by slightly touching the symbols on the b-control panel according to medical advice. Use a visual aid during application if needed.

Scope of application of b-intense® MED

Intended use

The **b-intense®** MED consisting of cabin shell, infrared back and front heater and accessories (speakers, ambient light) is intended for the local, contact free heat application in order to release muscular tensions and assist in the treatment of pain and tension conditions.

Application time: Depending on medical advice daily between 30 and 50 minutes.

Intended application

- › The **b-intense®** MED is intended for infrared therapy according to the specified purpose.
- › The product is exclusively intended for professional use under supervision of trained personnel and not for private use.
- › Intended application also includes fully reading and understanding the associated documentation and in particular these "safety instructions".

Not intended application

- › Every use that doesn't comply with the intended application is not intended and therefore prohibited.
- › Application without supervision enables undesired operating conditions during application and can harm the device or the user.
- › The use of not licensed components is prohibited. **b-intense®** is not liable for any damage occurring through unintended use. The risks of unintended use lies within the user.

Application specification

Intended medical indication:	Release muscular tensions and assist in the treatment of pain and tension conditions
Intended user population:	Intended user population are adult humans with an undisturbed sensitivity to heat. The intended user population do not have any contraindications mentioned in this IFU which excludes them from the heat treatment.
Intended part of the body or type of tissue applied to or interacted with:	Intended to interact with the back of the user
Intended user profile:	Physicians, registrated (RNs) nurses, Licensed practical (LPNs) Nurses, physical therapists
Intended conditions of use:	Stationary installation in a clinical or therapeutic environment.
Operating principle:	Local, contact free heat application with an electrical heater device on the back of the user

Operating conditions

Ambient Temperature: 18-35 °C

Relative humidity: 10-60 %

Atmospheric pressure: 70-106 kPa

General safety guidelines and use of the b-intense® MED

- › Read the Operation Manual thoroughly including the medical advice.
- › Don't use the **b-intense®** MED without supervision.
- › Always check the **b-intense®** MED for visual damage.
- › Ensure the conformity of the electric connection with the specification on the nameplate.
- › Ensure the proper connection and inspection of the shockproof socket,
- › Ensure the proper protection and connection to a fault-current circuit breaker (FI) of the socket.
- › Do not install the **b-intense®** MED in wet rooms, only in dry rooms (no condensing humidity).
- › Never cover the infrared heaters!
- › Keep easily flammable substances and fluids away from the infrared heaters.
- › Don't leave the **b-intense®** MED operating unsupervised.
- › In case of visual damage of the cabin contact **b-intense®**.
- › Never open the electrical components of the **b-intense®** MED.
- › In case of damage on the supply cable disconnect the **b-intense®** MED from the power supply and contact **b-intense®**.
- › Don't modify or adapt the **b-intense®** MED!
- › Exchange electronic components solely with original components of the manufacturer. All repairs must be performed by a qualified engineer authorized by **b-intense®**. The electrical components must not be opened.

Notes for safe use

- › Use the **b-intense®** MED according to the intended purpose defined in the operation manual.
- › The **b-intense®** MED must be used with naked, at the start dry back.
- › The **b-intense®** MED MUST NOT be used without headrest and backrest.
- › Open wounds, scars etc. must be covered.
- › The surrounding temperature should not be lower than 18 °C during application.
- › Preheat the **b-intense®** MED to 25–28 °C for ideal results.
- › Adjust the neck support so that the back of the head rests centrally.
- › Use less intensity (max. 50%) in the back heater in the beginning of the application for at least 5 minutes.
- › Raise the intensity of the back heater after heating up in small measures (5–10% per 5 minutes), so you always feel a pleasant warmth in the back.
- › Reduce the intensity of the back heater for the last 5–10 minutes of application about 10% per 2 minutes.
- › Enable enough air circulation between the back and the heater.
- › Be careful not to touch the surface of the heater with the back.
- › Watch out for a burning or painful feeling in the exposed area. Reduce the intensity if necessary.
- › The **b-intense®** MED is not intended for face application. Ensure to keep a safe distance to the front heater and do not look directly into the heater for more than 9 minutes. Use safety glasses if necessary (at least protection level 4-4 according to EN171).
- › Always use a towel for sitting.
- › Desinfect all areas which are in direct contact with the body according to the instructions after every application.
- › Clean the cabin regularly according to the instructions.
- › NEVER cover the infrared heaters or their ventilation system.

Pay attention to all safety notices and application notes for your own safety!

Warning

In case of any disturbances, damages or injuries prior to, during or after application, immediately report to b-intense®. For any questions or ambiguities please contact us at office@b-intense.at.

Medical mode of action

Infrared radiation (IR radiation) is part of the electromagnetic spectrum. The most important natural source of IR radiation is the sun. About 54 % of the solar radiation reaching the earth's surface is IR radiation. In addition, the human skin is increasingly exposed to artificial IR sources in the cosmetic and wellness sector, but also in the medical field. The application of electromagnetic radiation for heat application is widely spread.

Infrared rays are absorbed by the uppermost skin layers and converted into heat. IR exposure is therefore perceived as heat. Due to the thermal stress on the skin, the body reacts with heat defense reactions which are triggered gradually. The physiological response to heat depends on the tissue temperature, the duration and extent of the temperature increase and the size of the irradiated area.

The human organism is able to keep its core body temperature relatively constant regardless of fluctuations in ambient temperature. This is due to an efficient thermoregulation system, which acts as a negative feedback system to counteract deviations in the actual temperature value by more than ± 0.1 % of the setpoint. When heated, there is an increase in skin circulation and sweat secretion. As the vessels dilate, more blood can pass from the body core to the skin surface, the blood cools the skin and dissipates the heat. Deeper skin layers are not heated directly by infrared radiation, but by heat conduction. Because of subsequent evaporation at the skin surface the increased production of sweat has a cooling effect on the skin and is therefore an important thermoregulatory mechanism. In addition, sweat also contributes to the maintenance of optimal skin hydration.

Contraindications, warnings and precautions

Although every heat application bears the potential risk of damaging the skin, infrared radiation offers the advantage of contactless heat transfer, thus natural defense reactions remain largely unaffected.

Usually, unpleasant sensations or slight pain in the irradiated area, already clearly felt prior to burns, lead to behavioral changes (we move away from the source of radiation), so we are protected by our own natural averting reactions.

A coordinated response to heat stress of the entire cardiovascular system including skin circulation is necessary to respond appropriately to challenges to homeostasis. Users diagnosed with circulatory problems, cardiac conduction disorders including those with a pacemaker, malignant disease, epilepsy, and coagulation disorders have to be closely supervised and only treated after medical advice was sought.

Users with cancer and comorbid cardiac diseases have been found to be vulnerable to the short-term effects of extreme temperature.

Impairments in the thermoregulatory responses to heat stress are found in users with cardiovascular disease, such as congestive heart failure, or in individuals with type 2 diabetes mellitus.

The susceptibility to heat-related illness is also increased during use of certain medications (e. g. diuretics, antihistamines, antiepileptic drugs etc.), for users with Sjögren's syndrome and skin diseases (e. g. ichthyosis, anhidrotic ectodermal dysplasia), as in these cases sweating and thus the cooling effect on the skin is impaired.

Special caution is required with diseases that lead to increased metabolic activity and thus increased heat generation (status epilepticus, malignant hyperthermia).

Elderly people react particularly sensitively to higher temperatures. This also applies to children. The chance of overheating is particularly high for small children due to the large body surface area to body mass ratio.

In case of reduced or missing pain sensation, e. g. also under the influence of alcohol and drugs, tranquilizers and pain medication, thermal damage cannot be excluded.

In areas with extensive scars too little heat may be dissipated due to the reduced vascular supply. Special caution is also required with diseases that lead to increased metabolic activity and thus increased heat generation (status epilepticus, malignant hyperthermia). Therefore, scars in exposed areas should be covered.

Since the DNA repair mechanism may be impaired, we advise against infrared application for 24 hours after sunbathing or solarium. In case of acute injury/inflammation, the application of heat may lead to an increase in symptoms and should therefore be avoided.

In users with mild coagulation disorders, silent bleedings may occur, causing nevertheless inflammatory reactions, so heating is not recommended as possibly aggravating the condition.

Erythema ab Igne (EAI)

Repeated episodes of elevated skin temperatures can lead to pigmentary changes, which are known as Erythema ab Igne (EAI). Such reticulated, erythematous or hyperpigmented dermatoses have a good prognosis and mostly are a cosmetic problem and no serious damage. Reddened skin is normal to a certain degree and should decrease max. 2–3 hours after finishing the application. The use of the infrared cabin should be discontinued until skin has fully recovered. In case of persisting reddening of the skin (more than 12 hours), a doctor should be consulted. The skin structure could change permanently and heat melanosis could develop, should heat application nevertheless be continued. The use of the infrared cabin should be discontinued until skin has fully recovered.

The **b-intense**® MED is not intended for face application. A safe distance has to be kept to the front heater and directly looking into the heater for more than nine minutes is not recommended.

Information on radiation emissions

The **b-intense**® MED has one front heater and one back heater, which emit the following infrared radiation composition in normal operation and at maximum power:

- › Front heater: 16% IR-A ($0,78 \leq \lambda \leq 1,4 \mu\text{m}$), 51% IR-B ($1,4 \leq \lambda \leq 3 \mu\text{m}$) and 33% IR-C ($\lambda > 3 \mu\text{m}$)
- › Back heater: 15% IR-A ($0,78 \leq \lambda \leq 1,4 \mu\text{m}$), 58% IR-B ($1,4 \leq \lambda \leq 3 \mu\text{m}$) and 27% IR-C ($\lambda > 3 \mu\text{m}$)
- › The remainder is a minimal proportion of harmless visible radiation ($\leq 0,2\%$), ultraviolet or other proportions are not included.

Maximum output of OPTICAL RADIATION for all intended configurations of the **b-intense**® MED is 982 W/m^2
Spectral irradiance or spectral radiant exposure for all intended configurations of the **b-intense**® MED is 1000 W/m^2
Maximum variation of the output from the mean value across the treatment area is 50 W/m^2 .

By reducing the intensity of the infrared heaters, the wavelength components shift towards long-wave IR, which means that the components of IR-A and IR-B are less and IR-C more present. Provided that the body's natural reaction to move away from heat is not impaired, there is no possibility of burning the skin with the **b-intense**® heaters as long as they are used as intended.

Through suppression of heat sensibility (e. g. under the influence of drugs, alcohol or medication) or certain measures to increase blood flow in the skin (e. g. strong mechanical rubbing) burnt skin cannot be excluded. In such cases, medical supervision is necessary while using the cabin.

If used properly and as intended, damage to the eyes through IR-A and IR-B waves is very unlikely. The chance of damage can only grow in case of unfavorable gazing or sitting (gazing directly into the heater for more than 9 minutes or more often, especially with the upper body bent forward, so the eyes are relatively close to the heater). Should this not be ensured during application or should the body length of the user lessen the distance to the front heater, safety glasses are recommended (at least security level 4-4 according to EN 171).

The back heater may only be used with the provided backrest. Only then can the necessary distance to the heater and therefore the correct irradiance be guaranteed.

* Wavelength of electromagnetic radiation

Literature sources on medical effects are provided at: [b-intense.at/literatur-b-control-med](https://www.b-intense.at/literatur-b-control-med)

The current instructions (IFU) can be accessed via: **DE** <https://www.b-intense.at/ifu-med>, **EN** <https://www.b-intense.at/eu/ifu-med>

Previous versions can be inquired anytime via office@b-intense.at.

For questions or ambiguities please contact office@b-intense.at.


b-control MED: Operation

In the middle of the control panel is the digital display. The current values the cabin can be read in this field during the entire application as follows:

- › Heater Intensity
- › Session runtime in minutes
- › Current cabin interior temperature in degrees Celsius

The b-control MED unit

Starting up the b-intense® MED

When you touch the button  for the first time (approx. 1 second), the **b-intense®** MED switches on. By switching on, the following basic settings are automatically activated:

- › Runtime: The presetting of 30 minutes can be changed within the first 10 seconds after activation the **b-intense®** MED in 5-minute-steps using the +/- switch. The maximum limit is 50 minutes.
- › The light can be switched on or off using the light bulb button .
- › The back and frontheater activate in 2–3 seconds. The frontheater warms up to 100 % performance, the backheater to 50 % performance. The frontheater can be regulated via the operation panel. We recommend not to exceed 50 % in the backheater in the first 5 minutes to allow the body to adjust to the heat supply and the impact is increased.

	Function:	Operation:
	Backheater	Press the button to activate the intensity regulation for the backheater. Choose the desired intensity using the +/- buttons.
	Frontheater	Press the button to activate the intensity regulation for the frontheater. Choose the desired intensity using the +/- buttons.
	Light	Press the button to switch the light on or off.
	Music Playback	Press the button to activate the audio playback from the SD-card. Switch between the audio files using the +/- buttons.
	Volume	Press the button during audio playback to adjust the volume using the +/- buttons.
AUX	AUX-entry	Press the button once to activate audio playback from external sources, which can be connected via 3,5 mm audio cable to the AUX entry at the side of the power supply unit. This is the only feature of this button.

Note

All described functions on the control unit can be safely used by the user.


b-control MED: Service menu

The service menu allows adjustments to the basic settings after activation of the cabin. Switch the control off to arrive at the service menu. Press the following buttons in the given order.

1. 
2. 
3. 
4. **AUX**

After pressing the buttons in the correct order you get to the service menu.

- › Press the **Music button**  to switch between menu items.
- › Choose the distinct parameters with the +/- buttons.
- › Press the On/Off button to save the selection.

Parameters:

language: Choose between the standard languages English and German.

warm up: Adjust the warm up phase 0–50 minutes at 50% backheater performance. The default setting is 5 minutes.

start volume: The volume can be adjusted during application. The default volume is set at 30 of max. 100 and can be adjusted in the service menu.

runtime: The runtime can be adjusted for every individual application. Select the default runtime in the service menu (in steps of 5 minutes to max. 50 minutes).

automatic pre-heat: The cabin environment will be pre-heated automatically to the set pre-heat temperature after switch on. The treatment can be started at every time by pressing any button. The default setting is OFF.

pre-heat temperature: Adjust the cabin environment temperature for the pre-heat phase from 25°C to 40 °C. The default setting is 28 °C.

fan settings ON-AUTO-OFF: The fan function is not displayed on the control. The fan is always switched ON in default setting, switch it to OFF or AUTO if desired. The AUTO setting activates the fan automatically at 40 °C and deactivates at lower temperatures.

fan follow up time: After finishing the application, the fan will be activated for 10 minutes in case of higher temperatures than 40°C. This default setting can be adjusted to 5 or 15 minutes.

light follow up time: The diffused LED strip will still be activated for 10 minutes, so the cabin is not dark immediately after application. This can be adjusted to 5 or 15 minutes.

factory reset: The factory reset will set all parameters to default value

version: This panel holds information about the software version. This is used for traceability.

note:

Software updates must only be implemented by **b-intense®**.

Cleaning and care of the b-intense® MED

The **b-intense®** MED must be cleaned and disinfected after each use on parts in contact with the user.

Please note the warning notices for cleaning.

1. Unplug the infrared cabin from the main power supply before cleaning.
2. Use proper medical disinfectant for the surfaces. Mind the recommended exposure time.
3. Afterwards wipe it dry.
4. Ensure that no liquid enters the **b-intense®** heater.
5. The cabin must never be cleaned using a water hose, pressure washer, steam cleaner or splashing water.
6. The seats and backrests covered with imitation leather are easy to clean with a damp cloth. As the leather covers are made of special commercial artificial leather, they must be cleaned with a medical disinfectant for surfaces.
7. The protective grids are not to be cleaned.
8. To prevent the transmission of harmful germs, the door handles and the display must be disinfected after each use.

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Cleaning, disinfection and care

Regular maintenance of the **b-intense®** MED 110 V ensures proper and hygienic use. Pay attention to the following guidelines:

- › Position the **b-intense®** MED on a dry, even floor
- › Never clean the **b-intense®** MED wet.
- › Never use a pressure washer or a steam cleaner for external cleaning
- › If possible, use the recommended cleaning and disinfecting agents
- › Always use a dry towel underneath while seating

Recommended cleaning agents

- › Artificial leather: **b-intense®** recommends a mild, lukewarm soap sud.
- › Cabin walls: **b-intense®** recommends a mild, lukewarm soap sud.
- › Cabin ceiling and floor: **b-intense®** recommends a mild, lukewarm soap sud.
- › Glass: **b-intense®** recommends a standard glass cleaner without special features.
- › Metal components: **b-intense®** recommends a mild, lukewarm soap sud.

Use a lint-free microfiber cloth for cleaning. Never use abrasive or pungent detergents. Disconnect the **b-intense®** MED from the power supply for cleaning.

Control panel

Wet a microfiber cloth with a glass cleaner to clean the control panel. Dry the surface with a dry microfiber cloth. Never wet clean the control panel.

Infrared heaters

Before cleaning, ensure that the infrared heaters have cooled down completely. Wet a microfiber cloth with a glass cleaner to clean the infrared heaters. Dry the surface with a dry microfiber cloth. Never wet clean the infrared heaters.

Glass door

Use regular glass cleaner and a lint-free microfiber cloth for the glass components.

Disinfection

Disinfect all cabin parts, which regularly get in touch with the users body. This includes the upholstered parts, the control panel and the door handle (inside and outside) of the **b-intense**® MED. Use the surface disinfectant Lysoform Fugaten® Spray, which is medically approved and authorised by the manufacturer. Generously spray and leave it to dry according to the recommendations. After the exposure time, remove the remainder with a dry towel.

External cleaning

Clean the external areas of the **b-intense**® MED regularly. We recommend at least a weekly cleaning. Remove dust from the ceiling of the **b-intense**® MED at least monthly.

Inner cleaning

Clean the **b-intense**® MED regularly. We recommend at least weekly or daily, depending on the use. Use the recommended cleaners and cleaning aids. Additionally, use a soft brush for persistent contamination of the leather. The use of any abrasive products can lead to hard and brittle leather and is therefore discouraged.

Cleaning artificial leather

Wiping the material with a regular wet cotton cloth is sufficient for basic cleaning. In case of stronger contamination, a mild soap sud can be added, but should be tested on a hidden spot first. Remove the remainder with a cotton cloth and lukewarm water afterwards and dry the surface with a soft towel. Abstain from the use of products which contain grease, oil, dissolver, alcohol or abrasive products. Stains such as oil, grease, coffee, ink or textile discoloration have to be removed immediately, otherwise they can enter the material. The artificial leather should not be cleaned chemically. Only use fade resistant textiles on the artificial leather.

Note

Multiple cleaning and disinfection has no influence on the life time of the cabin part, only when used the recommended detergents and disinfectants. If not recommended detergents and disinfectants are used especially the surface of the artificial leather can break. Follow the instructions. If a visible defect is found on the cabin parts, please contact **b-intense**.

Maintenance instructions

The **b-intense**® MED requires no special maintenance. Annual preventive replacement (or replacement after 100 hours of operation) of two 10A fuses inside of the main unit by service personnel is recommended. The infrared lamps are not subject to aging and have full power until the end of their lifetime. If the upholstery of the bench, backrest or neckrest start to brittle through the years, they should be changed. Original spare parts can be ordered via **b-intense**®.

Disposal of the **b-intense**® MED

1. Make sure that the cabin has been disconnected from the power supply.
2. Dismantle the cabin in reverse order of the assembly instructions.
3. Dispose of the individual parts in accordance with the waste separation regulations (Waste Electrical and Electronic Equipment Directive WEEE Directive 2002 96/13/EC) in the respective country.
4. There are no extraordinary measures for decommissioning this device.

Certificate of Warranty

If unexpected faults or damage on the **b-intense**® MED 110V, you are covered for claims under the following guarantees:

- › The infrared lighting tubes come with a 10 year warranty starting on the date of purchase, provided that any defect is not due to external damage or handling by unqualified persons.
- › All other parts of the **b-intense**® MED come with a 2 year warranty, starting with the date of purchase, provided damage or faults are not related to external influences or handling by unqualified persons. Independent fitting or repairs, other than those specified in this manual, are excluded from the warranty. General wear and tear in principle is not covered by guarantee.
- › The scope of the guarantee for professionally installed replacement parts corresponds to that of the original warranty provided for the functioning of the cabin parts from the date of purchase. In addition, the minimum provisions apply under the statutory warranty from the date of installation.
- › Proof of purchase of the equipment must be provided with any claim under the guarantee by presenting the original receipt (delivery note or invoice).

b-intense® gmbh | Am Jungferenberg 17, A-2201 Gerasdorf by Vienna

These operating instructions have been subject to the official approval of an international verified inspection body. Information contained here is not a substitute for medical advice. The use of an infrared cabin for therapeutical purposes should only take place according to medical advice and under supervision from medically trained personnel.

Printing and typesetting errors excepted.

KEEP THESE INSTRUCTIONS SAFE FOR FUTURE REFERENCE AND HAND THEM OVER WITH THE B-INTENSE® MED IF THERE IS A CHANGE OF THE USER.

For further questions do not hesitate to contact us.
Your **b-intense**® team

b-intense®
INFRARED SAUNAS

b-intense® gmbh
branch office of CB-Capital
Beteiligungs GmbH
Am Jungferenberg 17
A-2201 Gerasdorf by Vienna, Austria

office@b-intense.at
www.b-intense.at

b-intense MED: marking plates

Main marking plate

Position of the label: outside, right side, near the mains connection

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SN 09/2020/2 P20202003-1
60 %
MD Class Ila
EN 60601-2-57:2011

Callouts:
- CE conformity mark
- Address of manufacturer
- Equipment type
- Mains supply data
- Serial number
- Follow instructions for use
- Separate collection electrical/electronic equipment

Marking plate front heater

Position of the label: back side of the front heater

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SN MEDBH 11-21-1234
60 %
EN 60601-2-57:2011

Callouts:
- CE conformity mark
- Address of manufacturer
- Equipment type
- Device connection data
- Serial number
- Follow instructions for use
- Separate collection electrical/electronic equipment

Marking plate back heater

Position of the label: back side of the back heater

The marking plate for the b-intense Back Heater contains the following information:

- Brand and Model:** b-intense® INFRAROTKABINEN
- Manufacturer:** b-intense gmbh, Zweigniederlassung der CB Capital Beteiligungs GmbH, Am Jungferenberg 17, A-2201 Gerasdorf bei Wien, Austria
- CE Marking:** CE 1304
- Equipment Type:** b-intense® Back Heater
- Technical Specifications:**
 - Ref. No. / Art. Nr.: 30100NM
 - Input voltage / Eingangsspannung: 230 V ~
 - Frequency / Frequenz: 50 Hz
 - Pmax: 300 W
 - Risk group / Risikogruppe: 2
- Warnings:**
 - Warning / Warnhinweis:** IR emitted from this device may cause eye injury. Avoid eye exposure. Aperture for optical radiation. / Die ausgestrahlte Infrarotstrahlung kann zu Augenverletzungen führen. Augenexposition vermeiden. Blende für optische Strahlung.
- Serial Number:** SN MEDBH 11-21-1234
- Additional Info:** 60%, EN 60601-2-57-2011
- Safety Symbols:** High voltage, Laser radiation, Warning, Grounding, Environmental protection, and disposal symbols.

Annotations on the image:

- Device connection data:** Points to the technical specifications section.
- Serial number:** Points to the SN MEDBH 11-21-1234.
- CE conformity mark:** Points to the CE 1304 marking.
- Address of manufacturer:** Points to the manufacturer's address.
- Equipment type:** Points to the b-intense® Back Heater text.
- Follow instructions for use:** Points to the warning symbols.
- Separate collection electrical/electronic equipment:** Points to the disposal symbols.

Marking plate b-control MED

Position of the label: metal cover mains unit

The marking plate for the b-intense b-control MED contains the following information:


- Brand and Model:** b-Intense® INFRAROTKABINEN
- Manufacturer:** b-intense gmbh, Zweigniederlassung der CB Capital Beteiligungs GmbH, Am Jungferenberg 17, A-2201 Gerasdorf bei Wien, Austria
- CE Marking:** CE 1304
- Equipment Type:** b-control MED
- Technical Specifications:**
 - Ref. No. / Art. Nr.: 40500
 - Input voltage / Eingangsspannung: 230 V ~
 - Frequency / Frequenz: 50 Hz
 - Pmax: 1000 W
- Serial Number:** SN 09/2020/1023
- Additional Info:** 60%, EN 60601-2-57-2011
- Safety Symbols:** High voltage, Laser radiation, Warning, Grounding, Environmental protection, and disposal symbols.

Annotations on the image:

- Device connection data:** Points to the technical specifications section.
- Serial number:** Points to the SN 09/2020/1023.
- CE conformity mark:** Points to the CE 1304 marking.
- Address of manufacturer:** Points to the manufacturer's address.
- Equipment type:** Points to the b-control MED text.
- Follow instructions for use:** Points to the warning symbols.
- Separate collection electrical/electronic equipment:** Points to the disposal symbols.

Marking plate packaging


Separate collection electrical/electronic equipment


The marking plate contains the following information and symbols:

- b-intense® INFRAROTKABINEN** (Manufacturer name)
- CE 1304** (CE conformity mark)
- b-intense gmbh, Zweigniederlassung der CB Capital Beteiligungs GmbH**
Am Jungferberg 17, A-2201 Gerasdorf bei Wien, Austria (Address of manufacturer)
- Type: b-intense MED** (Equipment type)
- Ref. No. / Art. Nr.: 12000**
- WEEE symbol** (Crossed-out bin): Separate collection electrical/electronic equipment
- Umbrella symbol**: Protection against moisture
- Glass symbol**: Fragile
- Upward arrows symbol**: This side up
- Pressure symbols**: 50 kPa, 110 kPa
- Temperature symbol**: 0° C, +40° C
- Humidity symbol**: 10 %, 100 %

Warning labels

b-intensive MED: Locations of EMISSION APERTURES

Hot surface. Do not touch the grid of the front and back heater.

If any malfunctions of the firmware are noticed, please unplug the b-intensive® MED and report immediately to office@b-intensive.

These operating instructions contain specific information about the radiation emissions, about options for radiation protection for users and about false use and risks through the installation (see chapter "Information on radiation emission").

b-intense MED: technical data

Mains voltage	230VAC
Frequency	50 Hz
Max. operating temperature	35° C
Safe work load	120 kg
Dimensions (W x D x H)	120 x 86 x 200 cm

Technical data b-control MED unit

Dimensions of power module:	300 x 210 x 60 mm
Dimensions of control unit (Front):	225 x 92 mm
Cut-out for the control unit:	215 x 82 mm (depth 20 mm)
Power module:	Protection class 1, Protection type IP20
Mains supply:	230VAC ±10% 50 Hz
IR Heater output:	2 x 230VAC max. 4A
Safety rating:	max. 16A, FI circuit breaker mandatory
Fan output:	12V max. 100 mA
Temperature reading:	2 x NTC
AUX IN:	Audio stereo input, 3.5 mm mini jack plugs
Speaker output:	max. 2 x 4 Watt RMS an 8 Ohm

Technical data LED light bar

Electrical data:

Operating voltage:	22–27V DC
Operating current:	max. 3 x 140 mA DC
Performance:	10 Watt
Protection class:	III connected to power supply according to EN 61347-2-13
Fuse:	0.25 A internal per colour channel
Reverse polarity protection:	up to 100
Connection cable:	LIYY 7 x 0.5 mm ² / 3.5 m

Case and dimensions:

Measurements:	Case (L x B x H) 600 x 100 x 70 mm
Shell:	aluminium white stove-enameled

Photometrical data:

Luminous flux (lm) white max.:	60 lm
Operating temperature max.:	50 °C

Technical data speakers

Size:	115 mm
Impedance:	3.4 Ohm ±10 %
Max. output:	40W

Connection and commissioning of the b-control MED

The following describes the installation and connecting steps of the b-control units. Please also note the attached assembly instructions.

1. Lay the cable for the front heater in the floor, the ends of the cables including their plug pins and the grey bow (= connection to the control unit FRONT/BACK) have to point to the control box and to the power module. These plugs will be plugged into the power module (control box). This handling ensures that the ends of the cables are connected correctly to the accurate socket of the front heater. (Please double check!) The cable of the front heater is marked with a red dot on the plug. This dot marker must be at the front!
2. The power module (control box) will be attached to the backside of the inspection board above the already prepared boreholes.
3. Plug the fan connection cable into the intended plug socket.
4. Plug the LED light bar into the intended plug socket.
5. Connect the control unit with the enclosed (RJ) cable as shown in the drawing.
6. The speakers are installed in the cabin ceiling and connected to the power module of the control unit via 2-pole plug. To enable the use of the audio input in the power module, the 3.5 mm jack socket on the power module must be connected with the enclosed cable. If the AUX input is not used, the AUX button has no function.
7. Connect the temperature sensor cable into the intended socket.
8. Connect the control unit with the enclosed (RJ) cable as shown in the drawing.
9. Plug the network cable into the power module.
10. Check the firm location of all connectors. Attention: To ensure sufficient ventilation, do not store other materials (insulation materials, instruction manuals, etc) inside the inspection flap.
11. Plug the power plug into a verified shockproof socket.

Possible error causes

The control unit cannot be switched on.

- no supply voltage, check house fuse
- mains plug unplugged
- Connecting cable from power module to the control unit is not plugged in properly

Control can be switched on, heaters do not warm up, temperature on the display >140 °C


- Temperature sensor defective or not plugged in

Control can be switched on, cabin temperature is displayed correctly but heaters do not warm up.

- Heater power is set to 0, adjust performance via control panel
- Heater is malfunctioning, contact **b-intense®**.

Should the error remain unresolved, please contact b-intense®.

Connection plan b-control MED control unit

Appendix A: Informations reg. EMF according to EN 60601-1-2

The **b-intense**® MED are intended for use in the electromagnetic environment specified below. The user of the above listed models should assure that they are used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The b-intense MED use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	

The **b-intense**® MED are intended for use in the electromagnetic environment specified below. The user of the **b-intense**® MED should assure that they are used in such an environment.

Enclosure Port			
Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±,8kV Contact ± 2,4,8,15 kV Air	±,8kV Contact ± 2,4,8,15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM fields and Proximity fields from RF wireless communications equipment IEC 61000-4-3	3V/m 80 MHz – 2,7 GHz 80 % AM 1 kHz	3V/m 80 MHz – 2,7 GHz	Mains power quality should be that of a professional healthcare facility environment.
	385 MHz (18Hz Pulse Modulation)	27 V/m	
	450 MHz (FM+/-5KHz deviation 1kHz sine or 18Hz Pulse Modulation)	28 V/m	
	710 MHz (217Hz PM)	9 V/m	
	745 MHz (217Hz PM)	9 V/m	
	780 MHz (217Hz PM)	9 V/m	
	810 MHz (18Hz PM)	28 V/m	
	870 MHz (18Hz PM)	28 V/m	
	930 MHz (18Hz PM)	28 V/m	
	1720 MHz (217Hz PM)	28 V/m	
	1845 MHz (217Hz PM)	28 V/m	
	1970 MHz (217Hz PM)	28 V/m	
2450 MHz (217Hz PM)	28 V/m		

	5240 MHz (217Hz PM)	9 V/m	
	5500 MHz (217Hz PM)	9 V/m	
	5785 MHz (217Hz PM)	9 V/m	
RATED power frequency magnetic fields IEC 61000-4-8	50 Hz or 60 Hz	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Comment: --			

The **b-intense**® MED are intended for use in the electromagnetic environment specified below. The user of the **b-intense**® MED should assure that they are used in such an environment.

Input a.c. power Port			
Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrical fast transient/bursts IEC 610004-4	± 2 kV 100kHz Repetition frequency	± 2 kV	Mains power quality should be that of a professional healthcare facility environment.
Surges IEC 61000-4-5	± 0,5kV, ±1 kV line(s) to line(s) ± 0,5kV, ±1 kV, ± 2 kV line(s) to ground (a)	± 1 kV Differential mode ± 2 kV Common mode	Mains power quality should be that of a professional healthcare facility environment.
Conducted RF induced by RF fields IEC 61000-4-6	3 Vrms 150 kHz - 80 MHz 6 Vrms ISM Bands c) 80% AM 1kHz	3 Vrms 150 kHz - 80MHz 6 Vrms ISM Bands c) 80% AM 1 kHz	Mains power quality should be that of a professional healthcare facility environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% UT; 0° 0% UT; 70 % 0% UT; 0%	0,5 Cycles 1Cycle 25/30 Cycles (50/60 Hz) 250/300 Cycles (50/60 Hz) (5s)	Mains power quality should be that of a professional healthcare facility environment. If the user of b-intense ® MED requires continued operation during power mains interruptions, it is recommended that the b-intense ® MED are powered from an uninterruptible power supply or battery.

Comment: a) Not applicable to CLASS II ME EQUIPMENT and ME SYSTEMS.

b) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

The **b-intense**® MED is intended for use in the electromagnetic environment specified below. The user of the **b-intense**® MED should assure that they are used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
Harmonic emissions IEC 61000-3-2	Class A	The model is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Appendix B: List of replaceable cables

Article Number	Description	Max. Length
80510	Power cable 230V	3 m

The logo for b-intense features the lowercase letter 'b' in a dark purple color. A dark grey, curved shape is positioned above the top of the 'b', resembling a stylized leaf or a wave. To the right of the 'b' is a registered trademark symbol (®).
b-intense®
INFRARED CABINS

www.b-intense.at